

The Specialty Excipient

Neusilin[®]

A totally synthetic magnesium aluminometasilicate (MAS) with exceptional excipient properties to improve API delivery and the quality of pharmaceutical preparations



Fuji Chemical Industries

Table of Contents

Characteristics	4
Chemical formula	4
General properties	5
Neusilin® grades	6
Typical properties	6
Typical applications and quantity needed in pharmaceutical preparations	7
Pharmaceutical applications	9
I. Oil and extracts to powder	10
II. Flowability Improvement	11
III. Anti caking	12
IV. Compressibility	15
V. Solid dispersion	16
VI. Hygroscopic velocity curve	16
VII. Hygroscopic equilibrium curve	17
VIII. Most compatible disintegrants	18
Cosmetic applications	18
Package size	19
Pharmacopoeia and regulatory information	
Dosage and safety	
Stability	
Conclusions	
Further reading	


Neusilin® - The specialty excipient

Neusilin® is a synthetic, amorphous form of magnesium aluminometasilicate. It is a multifunctional excipient that can be used in both direct compression and wet granulation of solid dosage forms. Neusilin® is widely used for improvement of the quality of tablets, powder, granules and capsules.


Neusilin® does not develop gel with aqueous solutions unlike other magnesium aluminum silicates.

The different grades of Neusilin® have been highly evaluated at home and abroad. It has a market presence of over 60 years in Japan.

Chemical formula




Neusilin® is amorphous and contains either tetrahedron or octahedron of Al, octahedron of Mg and tetrahedron of Si which are randomly attached to form a complex three dimensional structure. Neusilin® does not possess repeating units of a defined monomer.



Characteristics

1. Neusilin® is Magnesium Aluminometasilicate in either fine powder or granule form.
2. Neusilin® is represented by an empirical formula $\text{Al}_2\text{O}_3 \cdot \text{MgO} \cdot 1.7\text{SiO}_2 \cdot x\text{H}_2\text{O}$.
3. Neusilin® is amorphous, possesses very large specific surface area and has high oil and water adsorption capacity.
4. Neusilin® is superior in compressibility which enables to make hard tablets at low compression force. It can also improve hardness of other fillers and binders of low concentration.
5. Compounding with Neusilin® helps stabilize moisture sensitive as well as lipophilic APIs.
6. Neusilin® is stable against heat and has a long shelf life.
7. Neusilin® is available in various grades. The grades differ in their bulk density, water content, particle size and pH.
8. Neusilin® is an excellent carrier for solid dispersion via self-micro emulsifying drug delivery system (SMEDDS) and Hot-Melt Extrusion.

“Multi-problem solver”



General properties


Appearance	White powder or granules
Physical form	Amorphous
True specific gravity	2.0-2.2
Solubility	Practically insoluble in water and ethanol
Composition (%) on dried basis	Al ₂ O ₃ – 29.1 – 35.5 MgO – 11.4 – 14.0 SiO ₂ – 29.2 – 35.6
Loss on drying	Less than 20 to 5% depending on grades
CAS Number	12511-31-8
EINECS number	235-682-0

Neusilin® grades


UFL2	US2	S1	S2
Neutral	Neutral	Alkaline	Alkaline
Powder	Granule	Granule	Granule
Low water	Low water	High water	Low water

Electron micrographs


UFL2




US2



S1



S2



Typical properties

GRADE		Alkaline		Neutral	
		S1	S2	UFL2	US2
Appearance		White granule	White granule	White powder	White granule
Degree of whiteness (%)		>95	>95	>95	>95
Loss on drying (%) 110°C, 7 hours		13 - 20	< 5	<7	<7
Bulk density	Loose (g/ml)	0.30 - 0.37	0.29 - 0.37	0.06 - 0.11	0.13 - 0.18
	Tapped (g/ml)	0.36 - 0.43	0.34 - 0.42	0.10 - 0.17	0.16 - 0.22
True specific gravity		2.0	2.2	2.2	2.2
Specific surface area (m ² /g) ¹		110	110	300	300
Average particle size (µm)		112	115	3.1	106
Angle of repose (°)		30	30	45	30
Oil adsorbing capacity (ml/g) ²		1.3	1.4	2.7 - 3.4	2.7 - 3.4
Water adsorbing capacity (ml/g)		1.0	1.2	2.4 - 3.1	2.4 - 3.1
Acid consuming capacity (ml/g) ³		≥210	≥210	≥210	≥210
pH (4% slurry) ⁴		8.5 - 10.0	8.5 - 10.0	6.0 - 8.0	6.0 - 8.0

*1) BET surface area, nitrogen adsorption method

*2) Japanese Industrial Standard pigment test method (JIS K5101)

*3) Amount of 0.1N hydrochloric acid neutralized by 1g dried product (110°C, 7 hours)

*4) Weigh 2 g of sample, add water to make 50 ml. After stirring, allow to stand for 2 minutes, Measure pH using pH meter

Typical applications and quantity needed in pharmaceutical preparations

APPLICATION / FUNCTION	Quantity (%)			
	S1	S2	UFL2	US2
Diluent in solid dosage forms	30-90	30-90	30-90	30-90
Binder, increasing hardness, disintegration aid in tablets	5-20	5-20	1-10	1-10
Increase flowability	-	-	0.5-5	-
Anti-caking agent	-	-	0.5-5	-
Solidification of liquid API (eg: oil to powder)	-	-	30-50	30-50
For suspensions	-	-	1-5	-
Stabilization of deliquescent drugs	-	-	5-15	5-15
Solid dispersion, SMEDDS	-	-	20-50	20-50


Pharmaceutical applications

I. Oils and extracts to powder

Schematic flow



Oil adsorption capacity



Neusilin® US2 and UFL2 grades show higher oil adsorption capacity* when compared to MCC or colloidal silica.

*Linseed oil direct adsorption

Free flowing powder of linseed oil




Neusilin® US2 +30% linseed oil, Dry at 50°C



Linseed oil tablet, Ø11.3mm, 125N at 500 kg/cm²

Tablets of Scopolia extract and soybean oil


UFL2



A mixture containing 25% Scopolia extract or soybean oil and 25% UFL2 was compounded with equal amount of lactose. This mixture was subjected to static compression and tableting. We found no adhesion to pestle and mortar and the compressibility was good. The tablet did not exude the extract or oil during storage.

Tablets of Vitamin E


US2



An ethanol solution of tocopherol acetate (VE) 20-50% was compounded with proportional amount of Neusilin® and mixed well. To this mixture, 3% croscarmellose sodium and 1% magnesium stearate were added before tableting. High quality tablets with a load of up to 30% vitamin E can be prepared with Neusilin® US2.

II. Flowability Improvement



Angle of repose after adding Neusilin® and other excipients to potato starch





Neusilin® particles stick to the surface and aid flow



Potato starch only




Potato starch + Neusilin® UFL2



Electron micrograph showing Neusilin® UFL2 particles sticking to the starch surface. On addition to starch, the UFL2 particles stick to the surface and facilitate flow as in a 'roller blade' model. A 0.5% addition of UFL2 to potato starch vastly improves flowability.

III. Anti caking

Addition of 0.5% Neusilin® prevents caking UFL2



Neusilin® prevents caking at high humidity conditions UFL2



Sodium L-aspartate




Sodium L-aspartate + Neusilin® UFL2 0.5%

Condition: at 45°C, 75% RH, 2 days


IV. Compressibility

Neusilin® increases hardness of lactose tablet UFL2



Compounding lactose with 10% Neusilin® UFL2 results in higher hardness when compared to 15% microcrystalline cellulose

High quality tablets at low compression force US2



Tablet hardness of cornstarch/lactose based tablets compounded with either using Neusilin® US2, colloidal silica or calcium silicate. Corn starch, lactose and excipient were mixed thoroughly. Magnesium stearate as lubricant was added prior to tableting. Compression with Neusilin® US2 generally gives harder tablets compared to that with colloidal silica.

V. Solid dispersion

Formulating poorly water soluble drugs by solid dispersion leads to a remarkable improvement in dissolution and bioavailability. Neusilin® can potentially resolve problems associated with tableting and improve efficiency of solid dispersion.

Key advantages of Neusilin® as an adsorbent

- Flowability improvement
- High quality tablets at low compression forces
- High specific surface area
- High adsorption capacity
- Higher API load
- Restriction on reversion of amorphous form to crystalline state
- Inert core material

Case Study I


US2

Solid dispersion granules

Gupta MK, et al. Hydrogen bonding with adsorbent during storage governs drug dissolution from solid-dispersion granules. Pharm Res. 2002; 19:1663-72.



Dissolution profile of solid dispersion granules



Comparison of drug dissolution (after 30 min) from initial and stored solid-dispersion granules using USP Type II apparatus at 50 rpm. Data are shown for drug dissolution (% of initial) from solid-dispersion granules after storage at 40°C/75% RH (Gupta et al, 2002)

Case Study II

US2

Solid SEDDS (self-emulsifying drug delivery systems) formulation


Use of Neusilin® as adsorbent carrier to convert liquid SEDDS to solid SEDDS

1. Glyburide SEDDS tablets

Mura P, et al. New solid self-microemulsifying systems to enhance dissolution rate of poorly water soluble drugs. Pharm Dev Technol. 2012; 17:277-84

Self micro emulsifying formulation was prepared by adding under continuous stirring Tween 20 and Labrafac Hydro® WL (oil phase) and then distilled water to glyburide solubilized in Transcutol®. Glyburide tablets were prepared by direct compression. The preparation with Neusilin® US2 resulted in improved flow, compact tablets and improved dissolution profile.

Dissolution profile of glyburide preparation



Glyburide (GLY) dissolution profile from the different tablet formulations (Ref. Tablet – commercial GLY formulation; SME tablet – glyburide SME formulation consisting of Labrafac Hydro® as oil phase, Tween 20 as surfactant and Transcutol® as co-surfactant; TC tablet – glyburide formulation consisting of Transcutol® (TC) glyburide).


2. Solid SEDDS of paliperidone

Kanuganti S, et al. Paliperidone-Loaded Self-Emulsifying Drug Delivery Systems (SEDDS) for Improved Oral Delivery. J Disp Scie and Technol, 33:506–515, 2012

Optimized SEDDS formulation containing oleic acid, Tween 80 and Capmul® MCM L8 was adsorbed onto Neusilin® US2 to produce solid SEDDS (SEDDS-N). To understand the release behavior of paliperidone from solid SEDDS and pure drug, *in-vitro* dissolution test was performed.

Dissolution profile of paliperidone powder

The drug release was faster and the dissolution efficiency was higher for the solid SEDDS compared to that of crystalline form.



Case Study III US2

Hot melt extrusion (HME)


Maclean *et al.* Manufacturing and performance evaluation of a stable amorphous complex of an acidic drug molecule and Neusilin. *J Pharm Sci.* 2011; 100:3332-44.

HME of Sulindac-Neusilin® Drug Complex

Blends of Sulindac-Neusilin® in 1:1 and 1:2 (w/w ratio) were prepared by HME at 200°C.

Physical / Chemical Stability of Sulindac-Neusilin® HME Complex


The HME samples remained amorphous after 3 months of storage at 40°C/75% RH. The samples were found to remain amorphous for more than one year at ambient conditions.



Sulindac-Neusilin® HME tablets

Sulindac-Neusilin® 1:2 HME tablets showed 100% release in 90 minutes as against 9% release of Sulindac-Neusilin® crystalline tablets.


Dissolution profiles of HME Sulindac-Neusilin® tablets




VI. Hygroscopic velocity curve




1) 37°C, RH 53%




2) 37°C, RH 75%



3) 37°C, RH 92%




VII. Hygroscopic equilibrium curve



The hygroscopic equilibrium curve of different grades of Neusilin® indicate that Neusilin® absorbs very low amount of moisture up to 70% RH.


VIII. Most compatible disintegrants



The most compatible disintegrant with Neusilin® US2 was found to be croscarmellose sodium (Ac-Di-Sol) followed by cross-linked polyvinylpyrrolidone (Kollidon-CL) and carmellose calcium (ECG-505). The characteristics (large surface area and porous nature) of US2 and the cross linking of croscarmellose sodium act synergistically allowing the tablet to swell and absorb many times of its weight in water leading to quick disintegration. Neusilin® US2 improves flowability and makes sufficiently hard tablets at low compression forces. Increase in hardness and compression force did not affect the disintegration time or tablet conformity when croscarmellose sodium was used as a disintegrant.

As most of the starch-type disintegrants do not go well with Neusilin® US2, croscarmellose sodium is your best choice when you choose Neusilin® US2 in your formulations.

Cosmetic application UFL2




Unique properties of Neusilin® UFL2 make it an ideal component for cosmetic preparations. Although chemically the same as traditional crystalline Magnesium Aluminum Silicate (MAS), Neusilin® is both structurally and functionally very different. While MAS is used as a thickener, Neusilin® does not develop viscosity or form gel.

It is formulated in facial care products including lotions, eye shadow, cleansers, powders, acne and oily skin treatments and deodorants.



Neusilin® UFL2 applications



Oil adsorption

Neusilin® possesses excellent adsorption capacity for its extremely large specific surface area and porous reticulate structure. It can adsorb up to 330% of its own weight and maintains its powdery state. This superior adsorption capacity makes Neusilin® ideal for applications including oily and acne skin treatments.


Deodorization

Neusilin® works to eliminate odor in two ways. It can neutralize odor through making hydrogen bonds with bad odor compounds such as isovaleric acid or through covalent bonding with ammonia and trimethylamine. It can also physically absorb foul odor compounds by trapping them into its highly porous structure.

Functionalities

- Excellent sebum absorption capability
- High oil and water absorption
- Deodorant properties
- Pigment dispersion aid
- Anticaking agent
- Opacifying agent
- Ideal for powder lotions and other personal care products

Effective oil adsorption*



*Oleic acid adsorbed was washed with ethyl ether to determine the efficiency of adsorption.

Odor adsorption

Foul odor	Concentration / Amount used	Odor eliminated
Ammonia	1,000 ppm / 100 mg	79.5%
	1,000 ppm / 200 mg	96.2%
Trimethylamine	435 ppm / 100 mg	81.6%
Hydrogen sulfide	697 ppm / 100 mg	31.4%
Isovaleric acid	262 ppm / 20 mg	100%
Acetic acid	849.8 ppm / 100 mg	99.6%
Propionic Acid	382.3ppm / 20mg	100%
N-butyric Acid	311.8ppm / 20mg	100%

Package size

	Grade	Package size
Alkaline grades	S1	20 kg
	S2	20 kg
Neutral grades	UFL2	5 kg
	US2	10 kg

Samples are available upon request. Please contact your local distributor or sales representatives.

Pharmacopoeia and regulatory information

Neusilin® meets all the requirements of the current USP/NF and JPC. An US DMF type IV filed in 1998.

Dosage and safety

Neusilin® is safe with no reports of adverse reactions and is an accepted ingredient by the USP/NF and JPC. Based on the usage as an excipient in various formulations in Japan, Neusilin® up to 1.05 g can be used for oral uptake per day.* There are no established maximum oral intake limits specified by US-FDA.

*Encyclopedia of Pharmaceutical Additives, Japan, 2005

Stability

Neusilin® is a stable inorganic compound and meets JPC and NF specifications. The shelf life of Neusilin® is 3 years from the date of manufacture.

Conclusions

Neusilin® is a totally synthetic magnesium aluminometasilicate with exceptional excipient properties to improve API delivery and quality of oral solid-dosage form pharmaceuticals. Neusilin® is available in four grades and the two different pH options make it a versatile excipient for a wide-variety of applications.

With over 500 pharmaceutical preparations and a market presence of over 60 years in Japan, Neusilin® is well accepted by the formulators world-wide as an aid for formulations containing antibiotics, oily actives, poorly water soluble APIs, herbal mixtures, vitamins, etc. Neusilin® is also used as carrier for preparation of solid dispersion and self-micro emulsifying drug development systems.

Neusilin® has been demonstrated as an excellent adsorbent carrier for solid dispersion preparation via hot melt granulation, Self Micro-Emulsifying Drug Delivery Systems (SMEDDS) for BCS class II drugs such as meloxicam, naproxen, ketoprofen, glyburide and other highly permeable but poorly water soluble drugs. The most exciting use of Neusilin® is in Hot Melt Extrusion (HME). Neusilin® allows preparation of stable amorphous drug complex without any addition of polymers, waxes or plasticizers normally associated with HME. The samples can be recovered as amorphous powder and converted to highly stable tablets through direct compression.

Further reading

1. Shailesh T. *et al.*, Formulation and evaluation of liquisolid compacts for olmesartan medoxomil. *J Drug Deliv.* 2013;2013:870579.
2. Tan A, *et al.*, Transforming lipid-based oral drug delivery systems into solid dosage forms: an overview of solid carriers, physicochemical properties, and biopharmaceutical performance. *Pharm Res.* 2013;30(12):2993-3017.
3. Khanfar M, *et al.*, Mesoporous silica based macromolecules for dissolution enhancement of Irbesartan drug using pre-adjusted pH method. *Microporous and Mesoporous Materials.* 2013;173:22–28.
4. Janga KY, *et al.*, In situ absorption and relative bioavailability studies of zaleplon loaded self-nanoemulsifying powders. *J Microencapsul.* 2013;30(2):161-72
5. Kallakunta VR, *et al.*, Oral self emulsifying powder of lercanidipine hydrochloride:Formulation and evaluation. *Powder Technol.* 2012;221:375-382.
6. Sohn Y, *et al.*, Development of self-microemulsifying bilayer tablets for pH-independent fast release of candesartan cilexetil. *Pharmazie.* 2012;67:917-24.
7. Cha KH, *et al.*, Enhancement of the dissolution rate and bioavailability of fenofibrate by a melt-adsorption method using supercritical carbon dioxide. *Int J Nanomedicine.* 2012;7:5565-75.
8. Beg S, *et al.*, Development, optimization, and characterization of solid self-nanoemulsifying drug delivery systems of valsartan using porous carriers. *AAPS PharmSciTech.* 2012;13:1416-27.
9. Alam MA, *et al.*, Solid dispersions: a strategy for poorly aqueous soluble drugs and technology updates. *Expert Opin Drug Deliv.* 2012; 9:1419-40.
10. Kallakunta VR, *et al.*, Oral self emulsifying powder of lercanidipine hydrochloride: Formulation and evaluation. *Powder Technology.* 2012;221:375–82.
11. Kanuganti S, *et al.*, Paliperidone-loaded self-emulsifying drug delivery systems (SEDDS) for improved oral delivery. *J Disp Sci Tech.* 2012;33:506-15.
12. Hentzschel CM, *et al.*, Enhancement of griseofulvin release from liquisolid compacts. *Eur J Pharm Biopharm.* 2012;80:130-5.
13. Mura P, *et al.*, New solid self-microemulsifying systems to enhance dissolution rate of poorly water soluble drugs. *Pharm Dev Technol.* 2012;17:277-84.
14. Hentzschel CM, *et al.*, Suitability of various excipients as carrier and coating materials for liquisolid compacts. *Drug Dev Ind Pharm.* 2011;37:1200-7.
15. Hailu SA and Bogner RH. Complex effects of drug/silicate ratio, solid-state equivalent pH, and moisture on chemical stability of amorphous quinapril hydrochloride coground with silicates. *J Pharm Sci.* 2011;100(4):1503-15
16. Maclean J, *et al.*, Manufacturing and performance evaluation of a stable amorphous complex of an acidic drug molecule and Neusilin. *J Pharm Sci.* 2011;100:3332-44.
17. Boldhane SP and Kuchekar BS. Development and optimization of metoprolol succinate gastroretentive drug delivery system. *Acta Pharm.* 2010;60:415-25.
18. Chakraborty S, *et al.*, Utilization of adsorption technique in the development of oral delivery system of lipid based nanoparticles. *Colloids Surf B Biointerfaces.* 2010;81:563-9.
19. Hailu SA and Bogner RH. Solid-state surface acidity and pH-stability profiles of amorphous quinapril hydrochloride and silicate formulations. *J Pharm Sci.* 2010;99(6):2786-99.
20. Krupa A, *et al.*, Solid-state interaction of ibuprofen and Neusilin US2. *Thermochim Acta.* 2010;509:12-7.
21. Sander C, Holm P. Porous magnesium aluminometasilicate tablets as carrier of a cyclosporine self-emulsifying formulation. *AAPS PharmSciTech.* 2009;10:1388-95.
22. Vadher AH, *et al.*, Preparation and characterization of co-grinded mixtures of aceclofenac and Neusilin US2 for dissolution enhancement of aceclofenac. *AAPS PharmSciTech.* 2009;10:606-14.
23. Hailu SA and Bogner RH. Effect of the pH grade of silicates on chemical stability of coground amorphous quinapril hydrochloride and its stabilization using pH-modifiers. *J Pharm Sci.* 2009;98:3358-72.
24. Asai Y, *et al.*, Application of Neusilin UFL2 on manufacturing of tablets using direct compression method, Development of core tablets containing the function of small degree of decrease of hardness at the humid conditions. *Pharm Tech Japan.* 2009;25:67-70.
25. Asai Y, *et al.*, Application of Isomalt on manufacturing of tablets using direct compression method. *Pharm Tech Japan.* 2009;25:79-83.
26. Agarwal V, *et al.*, Dissolution and powder flow characterization of solid self-emulsified drug delivery system (SEDDS). *Int J Pharm.* 2009; 366:44-52.
27. Ruhland T, *et al.*, Nanoporous magnesium aluminometasilicate tablets for precise, controlled, and continuous dosing of chemical reagents and catalysts: applications in parallel solution-phase synthesis. *J Comb Chem.* 2007;9:301-5.
28. Bahl D and Bogner RH. Amorphization of indomethacin by co-grinding with Neusilin US2: amorphization kinetics, physical stability and mechanism. *Pharm Res.* 2006;23:2317-25.
29. Ito Y, *et al.*, Effect of adsorbents on the adsorption of lansoprazole with surfactant. *Int J Pharm.* 2005; 289:69-77.
30. Gupta MK, *et al.*, Formation of physically stable amorphous drugs by milling with Neusilin. *J Pharm Sci.* 2003; 92:536-51.
31. Gupta MK, *et al.*, Hydrogen bonding with adsorbent during storage governs drug dissolution from solid-dispersion granules. *Pharm Res.* 2002;19:1663-72.
32. Gupta MK, *et al.*, Mechanism for further enhancement in drug dissolution from solid-dispersion granules upon storage. *Pharm Dev Technol.* 2002;7:103-12.
33. Gupta MK, *et al.*, Enhanced drug dissolution and bulk properties of solid-dispersions granulated with a surface adsorbent. *Pharm Dev Technol.* 2001;6:563-72.

Neusilin[®]